KeyDenture Try-In™

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KeyDenture Try-In™

3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

Indications for Use: KeyPrint® KeyDenture Try-In™ is a liquid biocompatible photopolymer resin indicated for try-in dentures, intended to evaluate bite registration and occlusion for fully and partially edentulous patients.

Product Description: KeyPrint® KeyDenture Try-In™ is designed for additive manufacturing in vat polymerization 3D printers utilizing wavelengths between 385nm-405nm. KeyPrint KeyDenture Try-in™ is intended to be used within a computer-aided and manufacturing (CAD/CAM) digital dentistry system that includes a 3D scanner, design software, 3D printer, and post-cure unit. For any components that are used in conjunction with the KeyDenture Try-In™ resin, the user should review all applicable product labeling including Instructions for Use, user manuals, and other associated labeling. Strict adherence to all labeling is critical in assuring a safe and effective printed appliance.

Contraindications: Contains acrylate monomers and oligomers which, although rare, may cause an allergic reaction in individuals sensitive to acrylic containing products.

Warnings & Precautions:

1. Follow all compatible printer settings and curing times for safe and effective results. Refer to the QR codes.

2. Do not use any devices or components that are not KeyPrint® compatible or deemed acceptable by Keystone Industries. Refer to the QR codes below.

3. Review the product Safety Data Sheet (SDS) prior to use.

4. As per the SDS, wear proper personal protective equipment when handling KeyPrint® resins and uncured printed parts.

5. When pouring the resin, be careful not to splash.


6. Store in a cool, dry place 15°C-30°C (59°F-86°F) and away from light. Ensure that the bottle is capped while not in use.

7. In the unlikely event of a print failure during printing, filter the liquid resin through a mesh screen with pore sizing <200 microns. It is good practice to filter the resin in a vat periodically to prevent print failures.

8. Clean out the printer platform and vat prior to using a different batch of resin. Keystone recommends designating a vat that is specific for each KeyDenture Try-In™ shade. DO NOT mix different batches of the same product.

9. Keystone recommends against reclaiming the resin material without filtering.

Compatible Equipment: To ensure the biocompatibility of the final device, Keystone Industries collaborates with printer manufacturers to provide compatible printer and post-cure settings. Refer to the QR Codes for a full list of compatible printer and post cure units.

Compatible Printers

Compatible Post Cure Units

Directions for Use:

1. Ensure the resin is tempered to ambient temperature (20-25°C/68-77°F) prior to printing.

2. Agitate the bottle 1 hour prior to use for a minimum of 5 minutes to ensure all pigments are dispersed within the resin and to avoid bubbles during printing.

3. Once your design is completed per CAD software manufacturer’s directions for use, import the design into the CAM software unique to the printer manufacturer.

4. For optimal intaglio accuracy, avoid placing supports in the intaglio surface.

5. Printed parts should be cleaned with Isopropanol alcohol (> 97% purity) within approximately 8 hours from the completion of the print using the wash procedure below.

* Keystone discourages the use of denatured alcohol or ethanol for cleaning as they may diminish or degrade the quality of the finished parts.

Directions for post wash of printed part(s):

1. Using a soft brush or agitation, wash the printed parts in clean Isopropanol alcohol (>97% purity) bath for 30 seconds.

2. Remove printed parts from bath 1 and wash in another clean bath of Isopropanol alcohol (>97% purity) for an additional 30 seconds.

3. Remove printed parts from bath 2 and place in an ultrasonic bath for 2-3 minutes.

4. Remove printed parts from the ultrasonic bath and blow dry using compressed air.

5. Ensure all residual resin is removed from the parts, as residual resin will become shiny and highly noticeable. If resin is still present, use a soft brush to clean that isolated area.

Note: Printed parts should not exceed 5 minutes of IPA exposure as the parts may start to degrade.

Directions for post cure of printed part(s):

1. Refer to the QR codes for all compatible post cure settings.

2. For minimal surface defects or warping, allow parts to cool down before handling.

At this stage the medical device is cured and safe with respect to residual monomers.

Note: Prior to delivery to the patient, clean the medical device with soap and water to ensure the device is free of any debris from polishing or prior handling.

Disposal Considerations: Unused surplus and non-recyclable liquid resin materials may be considered hazardous waste, dispose of in accordance with all local, regional, national, and international regulations.

EN FOR THE FABRICATION OF TRY-IN DENTURES

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EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:


INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

keyprint®
by keystone Industries

KeyDenture Try-In™

keystoneindustries.com/symbolglossary

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

In case of Emergency:

INFOTRAC: 800-535-5053

Outside the US: (001) 800-535-5053

Customer# 76987

EN 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY

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