detax

3D Freeprint® Material

splint 2.0

Medical resin for hard splints

bg	Инструкция за употреба	5
cs	Návod k použití	8
da	Brugsanvisning	11
de	Gebruiksaanweisung	14
el	Οδηγίες χρήσης	17
en	Instructions for use	20
es	Modo de empleo	23
et	Kasutusjuhend	26
fi	Käyttöohje	29
fr	Mode d'emploi	32
hr	Upute za uporabu	35
hu	Használati útmutató	38
it	Istruzioni per l'uso	41
lt	Naudojimo instrukcija	44
lv	Lietošanas instrukcija	47
nl	Gebruiksaanwijzing	50
no	Bruksanvisning	53
pl	Instrukcja stosowania	56
pt	Manual de instruções	59
ro	Instrucţiuni de utilizare	62
ru	Инструкция по применению	65
sk	Návod na použitie	68
sl	Navodila za uporabo	71
sv	Bruksanvisning	74
tr	Kullanım kılavuzu	77
Annex 1		80

Ordering information

Freeprint® splint 2.0 385 clear

500 g	02080
1000 g	02076

Freeprint® ibt 385 clear

1000 g	04249
--------	-------

Freeprint® crown 385 500 g

A1	02372
A2	02378
A3	02417
A3,5	04709
B1	02481
B3	02645
C2	02766
D3	02783
BL	02845

Freeprint® crown 385 1000 g

A1	02376
A2	02415
A3	02446
A3,5	04708
B1	02519
B3	02758
C2	02782
D3	02825
BL	02884

Ordering information

Freeprint® denture impact 1000 g

dark reddish pink	04727
pink	04437
pink-transparent	04436

Freeprint® tray 2.0 385 green

1000 g	02505
5 kg Ecobag	04624

Freeprint® ortho 385 clear

1000 g	03989
5 kg Ecobag	04323

Freeprint® ortho 405 clear

1000 g	03988
--------	-------

Freeprint® temp 385 500 g

A1	04058
A2	04059
A3	04060

Freeprint® temp 385 1000 g

A1	04062
A2	04063
A3	04064

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Resin for dental 3D printing

Indication

Splints, fixation and transfer keys

Patient Target Group

Persons being treated in the context of a dental procedure.

Intended Users

Dentist, dental technician

Suitable for the following DLP printers/cleaning/post-exposure

see "Annex 1"

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Minimum material thickness for the design is 1.5 mm
- Polish surface mechanically.
- Processing temperature (at 23 °C ± 2 °C / 73 °F ± 4 °F).

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.
- Biocompatibility is only guaranteed with complete polymerisation.
- MMA-free.
- Refer to the relevant safety data sheet for hazard and safety information.

Cleaning

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards.

Note

- Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- Do not use any acids or solvents.

Notes

- Detax shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Please follow the instructions on the safety data sheet!

For users and/or patients

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- Store Freeprint® splint 2.0 in a dry place 15 °C - 28 °C / 59 °F - 82 °F and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of Freeprint® splint 2.0 may cause allergic reactions in predisposed persons. In such cases refrain from using the product. Freeprint® splint 2.0 only insert intraorally in completely polymerised state.


Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage

Application

At 23 °C ± 2 °C / 73 °F ± 4 °F

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

① Construction process

Generation of a Print Job complying with machine and material parameters

② Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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*The set of parameters includes all relevant material- and printer specific information

Annex 1

Digital light Processing (DLP) printer, operation software and parameter

Model (Picture)	Printer Manufactur- er, Model	Light source	Light intensity	Operation Software	Parameter data set*
	Microlay Versus	385 nm	4.3 mW/cm ²	Microform 1.0.3.7	DETAX Freeprint Splint 2.0 385100 microns v5.3
	Miicraft Alpha	385 nm	6.0 mW/cm ²	Utility 6.4.4	DETAX freeprint splint 2.0
	Miicraft Prime Series	385 nm	6.2 mW/cm ²	Utility 6.3.0	Detax_Freeprint splint 2.0 385_100_PH
	Miicraft Ultra Series	385 nm	5.7 mW/cm ²	Utility 6.3.0	Detax_Freeprint splint 2.0 UV_100
	Rapidshape D20+/D30+/ D40+	385 nm	2.0 mW/cm ²	Netfabb 2020	DETAX Freeprint splint 2.0 385
	Way2 Production SolFlex Series	385 nm	8.0 mW/cm ²	Netfabb 2020	Freeprint splint 2.0 385

*The set of parameters includes all relevant material- and printer specific information

Cleaning Equipment

Cleaning unit

Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS Wash	Use the following settings: DETAX Freeprint-splint 2.0 uv. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

Light curing Equipment

Light Curing unit

Manufacturer, Model	Curing process
---------------------	----------------

Dentalfarm Photopol	2 × 3 min, progressive + N ₂ , turn around components after 3 min
------------------------	--

Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 100 %; Duration = 120 s Post curing B: Wavelength = 460 nm; Intensity = 100 %; Duration = 120 s
-----------------------------	--

NK Optik Otoflash G171	2 × 2000 flashes under inert gas, turn around components after 2000 flashes
---------------------------	---

NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz
------------------------------	-------------------------------------

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

detax GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen · Germany
+49 7243 510 0 · post@detax.com · detax.com

07/2026
3